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Keto-Diastix® Reagent Strips

Test for Glucose and Ketone (Acetoacetic Acid) in urine.

SUMMARY AND EXPLANATION/INTENDED USE: KETO-DIASTIX® Reagent Strips are intended for self-testing by persons with diabetes and for healthcare professionals to monitor for the presence and concentration of glucose and ketone (acetoacetic acid) in urine.^{1,2,3} KETO-DIASTIX Strips are specific for testing for glucose and ketone in urine. Acetoacetic acid, more commonly referred to as "ketone body," can be found in the urine from persons with diabetes. When dipped in urine, these test areas on the KETO-DIASTIX Strip change color according to the amount of glucose and ketone in the urine. Use of KETO-DIASTIX Reagent Strips can alert you and your doctor to changes in your condition for which adjustments in your diet and/or medication may be needed. Carefully follow the testing schedule your doctor or diabetes educator establishes for you. Normally, the body burns glucose (sugar) for energy; however, because of the absence of or an insufficient amount of insulin in people with diabetes, glucose cannot be burned to supply energy for the body. Much of the glucose is expelled in urine. Your body turns to the fat stored in your body to get its energy. When this happens, ketones are formed in the blood and finally go into the urine. The presence of ketone is really the body's signal that something out of the ordinary is happening. The presence of glucose and/or ketone in the urine can be an indication your diabetes may be out of control. Small amounts of ketone in your body are usually not harmful. However, when you have too much of this substance in your body, it can be dangerous. If you burn too much fat and produce a lot of ketones, you may develop a complication called diabetic KETOACIDOSIS. This can lead to diabetic coma. You need to take fast action when you see the ketone warning signal. You should contact your doctor or diabetes educator immediately. Your doctor or educator can tell you how to change your insulin dose or your food intake to bring your diabetes back into control and to prevent ketoacidosis and coma.

Read the KETO-DIASTIX insert carefully and completely before you begin testing.

WHEN SHOULD YOU TEST WITH KETO-DIASTIX REAGENT STRIPS?

1. When you have a cold, the flu or any other kind of illness.
2. When you "feel" the signs of high blood sugar (over 13.3 mmol/L or 240 mg/dL) or when your blood sugar is well over the range your doctor has set for you (if you do blood glucose monitoring).
3. When you are under stress. (e.g., at home, at school, at work, or on a vacation).
4. Regularly during pregnancy or after a testing pattern has been established with your doctor or educator.

SPECIMEN COLLECTION AND PREPARATION: Collect urine in a clean, dry container and test it as soon as possible, or alternatively pass the Reagent Strip through the urine stream. If testing cannot be done within an hour after collecting a urine specimen, tightly cover the specimen and refrigerate it immediately. Let it return to room temperature before testing. Prolonged exposure of unpreserved urine to room temperature may result in microbial contamination and bacterial consumption of urine glucose. This may affect the accuracy of the result.

Urine preservatives may affect test results. Some preservatives do not adequately protect glucose from being metabolized by contaminating or infecting organisms.

STORAGE AND HANDLING: Never use the Reagent Strips past the expiration date printed on the bottle label. A new bottle of test strips can be used for 6 months after being first opened. Always write the date you first opened the bottle on the bottle label. Do not use product (opened or unopened) after expiration date. Use of strips beyond the expiration date may yield low results.

To keep KETO-DIASTIX Reagent Strips fresh, follow these suggestions carefully.

- Keep your finger or other objects from touching the reagent areas before testing.
- Store at room temperatures between 15°–30°C (59°–86°F).

- Do not store the bottle in direct sunlight.

IMPORTANT: PROTECTION AGAINST AMBIENT MOISTURE, LIGHT AND HEAT IS ESSENTIAL TO GUARD AGAINST ALTERED REAGENT REACTIVITY.

- Discoloration or darkening of reagent areas may indicate deterioration. If either of the reagent areas are ever discolored or darkened, throw the strip away and use a strip from a new bottle.
- Keep unused Reagent Strips in the original bottle with the cap tightly closed.
- Always replace the cap immediately and tightly after removing a Reagent Strip.
- Never transfer Reagent Strips to another bottle.
- Do not remove desiccant from bottle. The desiccant absorbs moisture and keeps the strips dry.

WARNINGS AND PRECAUTIONS: KETO-DIASTIX Reagent Strips are for *in vitro* [iv] diagnostic use only. *Do Not swallow.*

DIRECTIONS FOR TESTING:

Materials needed:

- KETO-DIASTIX® Reagent Strips
- Timer or watch (capable of measuring in seconds)
- A clean, dry container

Step 1. Collect a fresh urine specimen in a clean, dry container. (Or you may choose to pass the test end of the strip through a stream of urine.)

Step 2. Open the bottle and remove one Reagent Strip. Immediately replace cap tightly on bottle. Hold the plastic end of the strip. DO NOT touch the test areas on the strip.

Step 3. Dip the test areas of the strip into the urine and remove immediately, drawing the edge of strip against the rim of the urine container to remove excess urine. Or you may pass the test end of the strip through a stream of urine. Begin timing.

Step 4. Exactly 15 seconds after dipping compare the ketone test area to the Ketone Color Chart.

Step 5. Exactly 30 seconds after initial dipping compare the glucose test area to the Glucose Color Chart.

NOTE: The total time for glucose and ketone tests together is 30 seconds. Ignore any color changes that may occur after the specified reading times.

Step 6. Record results of both tests.

IMPORTANT: If the test shows a moderate or large amount of ketone, call your doctor or diabetes educator.

EXPECTED VALUES:

Glucose: Small amounts of glucose are normally excreted by the kidney.⁴ These amounts are usually below the sensitivity of this test but on occasion may produce a color between the negative and the 5,5 mmol/L (100 mg/dL) color blocks. Results of 5,5 mmol/L (100 mg/dL) may be significantly abnormal if found consistently.

Ketone: Normal urine specimens ordinarily yield negative results with this test. Detectable levels of ketone may occur in urine during physiological stress conditions such as fasting, pregnancy and frequent strenuous exercise.^{2,3} In ketoacidosis, starvation, or other abnormalities of carbohydrate or lipid metabolism, ketones may appear in urine in large amounts before serum ketone is elevated.⁴

RESULTS: Results with KETO-DIASTIX Reagent Strips are obtained directly from comparison to the Color Chart. The color blocks represent nominal values; actual values will vary around the nominal values.

Vertrieb durch:
Distribué par :
Gedistribueerd door:
Ascensia Diabetes Care NV-SA
Pegasuslaan 5
BE-1831 Diegem (Machelen)
Tel. 0800 97 212 (BE)
Tel. 800 27 254 (LUX)
info@ascensiadibetescare.be
www.ascensiadibetescare.be
www.diabetes.ascensia.com

Vertrieb durch:
Ascensia Diabetes Care Deutschland GmbH
51355 Leverkusen
Deutschland
Ascensia Diabetes Service
Telefon: 0800 7261880 (kostenfrei)
E-Mail: info@ascensia.de
www.diabetes.ascensia.de
www.diabetes.ascensia.com

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Ascensia Diabetes Care Holdings AG
Peter Merian-Strasse 90
4052 Basel, Switzerland

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© 2019 Ascensia Diabetes Care Holdings AG
Ascensia Diabetes Care Australia Pty Ltd
Ascensia Diabetes Care Hong Kong Limited